



Manufactured
in the UK



**Safelab Systems
Products**

Filtration Fume
Cupboards

Ducted Fume
Cupboards

Extract Systems

Cyclonic Scrubbers

Horizontal Laminar
Flow Cabinets

Vertical Laminar
Flow Cabinets

Safety Cabinets

Forensic Safety Cabinets

Polypropylene Cabinets

Glove Boxes

Storage Cabinets

**Safelab Systems
Service**

Maintenance Contracts

Airflow Testing

Ductwork Testing

Filter Testing

DOP Testing

KI Discus Testing

Installation

Commissioning

Replacement filter
service

Airone 1000PPCSC Chemical Storage Cabinet



Designed specifically for the safe storage of chemicals.

Protects staff and the environment.

Benefits:

Quiet operation/running mode

Removable perforated shelves

Large capacity carbon filter for the removal of Formaldehyde

Optional internal lighting available

Dimensions:

External (WxDxH) 1000 x 700 x 2140mm

Internal (WxDxH) 980 x 675 x 1530mm

Also available in 1200 width (external)

Construction & control system:

White polypropylene body & shelves

Clear PVC door

Mains On/Off

Fan failure alarm

Safelab Systems Ltd

Unit 29 Lynx Crescent Weston Super Mare North Somerset BS24 9BP

Tel: +44 (0) 1934 421 340

Fax: +44 (0) 1934 641 569

E-mail: safelab@safelab.co.uk

www.safelab.co.uk